

# Programmatically Determine DAQmx Device Name when Deploying Executables

---

## Introduction

When preparing a program to be built into an executable meant for deployment, it may be advantageous to include code that will detect all devices on the deployment machine, determine which one is the correct device, and then output its DAQmx resource name for use within the executable. This will cut down on confusion for any operators who may be running the executable and unfamiliar with LabVIEW or the DAQmx device resource naming paradigm.

The example program downloadable here uses a DAQ System Property node to return an array containing the DAQmx resource names of all of the hardware on the computer. It then uses a for loop to iterate through each resource name, obtaining the Product Type and Serial number of that specific device. These can be compared to the desired value for the device that was deployed with the code. When a match is found, the for loop is exited and the DAQmx Device name is output for use in the acquisition code.

Note: If you have multiple devices of the same type on the computer, this code will stop after it comes across the first one. In the case of multiple devices, it may be best to use the serial number as the identifier instead of the product type.

## Steps to Complete

1. Enter the Product Type of the device you want to acquire from, ex USB-6009 or PXI-4472
2. Enter the channels of the device you want to acquire from, beginning with a "/", ex "/ai0:3" to acquire from analog input channels 0, 1, 2, and 3.
3. Configure the Maximum and Minimum expected values of the acquisition, as well as the desired Timing Properties.

## Programmatically Determine DAQmx Device Name when Deploying Executables

4. Run the VI. It will now search the system for the matching product and acquire on the configured channel.

For instructions, select **Help>>Show Context Help**

ProductType you are looking for  
PXI-4472 e.g. USB-6009

DAQmx Device Name  
1/0 Device Name of the desired Device  
e.g. Dev1

Desired Channels  
/ai0:3 These are the channels that you want to acquire on. These are concatenated with the device name and must therefore start with a "/"  
e.g. /ai0:3

**Channel Parameters**

Minimum Value -10.00 Physical Channels This is an indicator. It is the result of the name configured in the box above

Maximum Value 10.00

**Timing Parameters**

Samples per Channel 1000

Sample Rate (Hz) 10000.00

Measurement

Amplitude

10.0  
8.0  
6.0  
4.0  
2.0  
0.0  
-2.0  
-4.0  
-6.0  
-8.0  
-10.0

10:37:53.96 PM

DevNames

0

All device names present in the system, should look like MAX's DAQmx Devices list



